



Cleaaver

S C I E N T I F I C



Fermentation and Bioprocess System



Benchtop
0.5 to 20 L



Pilot Scale
30L & 50L



Production Scale
>100L



Topics



► Benchtop fermentation system

- Controller
- Vessel Types
- Optional Devices & Accessories
- Fermenter Interface & Features

The 4th Generation Controllers

FS-05



FS-06




FS-07



▶ FS-05 Parallel Fermentation System



- 10.4" Touch Screen
- Compatible Vessel Volume: 0.5L, 1L, 3L, 5L, 10L, 15L, 20L
- Connect up to two vessels and compatible with all types of vessels
- Standard METTLER TOLEDO pH & DO sensors
- Eight built-in peristaltic pumps
- Optional devices available & 4 external peristaltic pumps



▶ FS-06 One Fermentation System



- 8" Touch Screen
- Compatible Vessel Volume: 0.5L, 1L, 3L, 5L, 10L
- Three peristaltic pumps
- Standard **METTLER TOLEDO** pH & DO sensors
- Only **1** external peristaltic pump

#COMPACT
#SPACESAVING
#PERSONALUSE



► FS-07 Evo Fermentation System



- 10.4" Touch Screen
- Compatible Vessel Volume: 0.5L, 1L, 3L, 5L, 10L, 15L, 20L
- Compatible with all types of vessels
- Standard **METTLER TOLEDO** pH & DO sensors (Dissolved Oxygen)
- Four built-in peristaltic pumps
- Optional device available & 2 external peristaltic pumps

#COMPATIBILITY



▶ **FS-06-EPM**





▶ **Expansion Module for FS-06**



- Flexible and expandable
- Additional Software Feature: 2-stage DO cascade, pH & DO Stat
- Optional Device Available: 2 external peristaltic pump
 - Oxygen Enrichment Module
 - Gas Mixing Station
 - Photosynthesis Lighting Module
 - ORP Probe

► Controller Selection Guide



				
Model	FS-05	FS-06	FS-06-EPM	FS-07
I want to add additional pumps to the controller	✓	✓	✓	✓
I want to work with small pilot-scale, 15L or 20L	✓			✓
I have limited lab space		✓	✓	
I want to do a parallel study	✓			

► Specifications



Model	FS-05	FS-06	FS-06 + FS-06EPM*	FS-07
Heating System	Duo heating			
Working Volume Range	500ml - 20L	500ml - 10L	500ml - 10L	500ml - 20L
Autoclavable Glass Vessels	Yes			
Interchangeable Vessels	Compatible with all types of vessels, except for 5L solid state which is only for FS-05 and FS-07 controller			
Number Of Vessels Controlled Per Controller	2	1	1	1
Number Of Vessels Controlled Via Remote Software	Max 32	Max 16	Max 16	Max 16
Touchscreen Controller	10.4"	8"	8"	10.4"
Number Of Peristaltic Pumps	8	3	3	4
Gas Mixing Options	Available	No	Available, *	Available
Oxygen Enrichment	Available	No	Available, *	Available
Mass Flow Controller	Available	No	No	Available
Off Gas Analyzer	Available	No	No	Available
ORP Probe	Available	No	Available, *	Available
Lighting Module	Available	No	Available, *	Available
External Pump	4 max.	1 max.	2 max.	2 max.
Solid State	Available	No	No	Available

* Optional expansion module (FS-06-EPM) needed.




► Features

- Five interchangeable types of glass vessel
- Submerged and solid state integration
- Intuitive user interface
- 15-step sequence program
- Multilanguage user interface
- EZScript for advanced process control
- Smart communication technology

► Features



► Five Interchangeable types of glass vessels



*Compatible with all types of vessel, except 5L solid state (only usable with FS-05 and FS-07)



► Features

➤ Submerged and solid state integration



**Submerged
Culture**



Solid State



**Submerged
Culture**



*Compatible with all types of vessel, except 5L solid state (only usable with FS-05 and FS-07)


► Features

➤ Variable-speed Peristaltic Pumps




- Dual direction 
- Variable speed: 1-65 rpm
- Manual or 15-step sequence control

Fed-batch fermentation





Continuous culture



► Features

- Intuitive user interface
- Multilanguage user interface



► Features

➤ 15-step sequence program

- Choose either **Manual** or **15-step sequence control** for Agitation, Temperature, O₂, Feeding, etc.

Agit Program

	1	2	3	4	5	6	7	8
RPM	100	150	200	300	100	0	300	600
MIN.	60	50	80	0	120	0	50	30
Now	➡	⬅	⬅	⬅	⬅	⬅	⬅	⬅
	9	10	11	12	13	14	15	
RPM	800	500	500	300	200	150	300	
MIN.	200	60	60	120	60	10	300	
Now	⬅	⬅	⬅	⬅	⬅	⬅	⬅	

Prev. Step Next Step Close



Feed1 Program

	1	2	3	4	5	6	7	8
RPM:	20	45	60	30	65	45	40	0
ON:	60	2	5	5	10	15	8	0
OFF:	60	0	2	2	5	6	6	0
Cycle:	360	30	5	0	100	85	16	0
Now	➡	⬅	⬅	⬅	⬅	⬅	⬅	⬅

	9	10	11	12	13	14	15
RPM:	0	0	0	0	0	0	0
ON:	0	0	0	0	0	0	0
OFF:	0	0	0	0	0	0	0
Cycle:	0	0	0	0	0	0	0
Now	⬅	⬅	⬅	⬅	⬅	⬅	⬅

ON/OFF unit: second

Prev. Step Next Step Close

► Features

➤ EZScript for advanced process control

- Write your own program using a computer language called BASIC
- Almost everything can be controlled



Overview	Control	Pumps	Calibration	Charting	System Setup	Service
Scripts			Console			
Print "your batch details"			Your batch details			
Agit_SV (300)			Set Agitation 300rpm			
Temp_SV(27.8)			Set Temp 27.8C			
PH_SV(5.91)			Set pH 5.91			
DO_SV(40)			Set DO 40%			
If(DO_PV (<40) MANU_DO_CASCADE agit Agit_DOCA_min (100) Agit_DOCA_max (500) Agit_DOCA_step (10) Check DO_PV (= 40)			If DO is lower than 40% Increase agitation speed From 100rpm To 500rpm With 10rpm increment Check if DO is 40%			


► Features

➤ Smart communication technology

Control up to 16 systems via remote monitoring



#COMPLIMENTARY




► Features

➤ Smart communication technology



1 PC
16 Controllers
32 Vessels



Topics



- ▶ Benchtop fermentation system
 - Controller
 - Vessel Types
 - Optional Devices & Accessories
 - Fermenter Interface & Features

► Vessel Types



Type A

Jacketed vessel



Type B

Single wall



Type B
Blanket

B+ blanket. Single wall with blanket



Type C

Air lifter



Type D


Single wall with heating base

Solid state



Solid state

► Typical Stirring Tank Design




Medium-free: Sanitary SUS 316

**Premium German-made
Borosilicate Glass**

**2:1 or 3:1 (air lifting) H/D Ratio
ideal for
mixing & homogeneity**

**In contact with medium:
Sanitary SUS 316**

► Typical Stirring Tank Design



Impeller: Mixing & Breaking up bubbles



Rushton
Standard




Pitched Blade
Shear-force
sensitive cells



Customized
impellers
available

► Typical Stirring Tank Design

Headplate arrangement




Customization
upon request

▶ FS-V-A Type Vessel-**Thermostat** Double Jacketed



▶ Key features:

- Uniform temperature control
- Suitable for temperature and **shear-force**  **sensitive cells** such as **mammalian** and **insect cell cultures** (with pitched blade impellers)



FS-V-AS5L




FS-V-A01L

► FS-V-B Type Vessels – **Thermostat** Single Wall Dish Bottom



► Key features:

- Stainless steel cooling coil for temp. control
- No dead volume with dish bottom
- Suitable for **microbes**



▶ FS-V-B Type Vessels – Dry Heating Single Wall w/ Heating Blanket




▶ Key features:

- Heat evenly and quickly
- Ideal for **photosensitive** and **photo-inhibition** cultures, such as **microalgae** and **nitrifying bacteria**
- Suitable for cultures which can withstand quick temp. change



FS-V-B05L

► FS-V-B Type Vessels – Dry Heating Single Wall w/ Heating Blanket



FS-V-B05L

▶ FS-V-C Type Vessels –Air Lifter Thermostat & no temperature control

▶ Key features:

- No impeller design: micro-sparger and adjustable draft tube to facilitate water circulation
- Single or Double -Walled 3:1 glass vessel
- Suitable for shear-sensitive cells and photosynthesis reaction: **plant cell, microalgae, cyano-bacteria**



FS-V-C053 FS-V-C054




Lighting
Module
(FS-O-PB-2)

► FS-V-D Type Vessels – Dry Heating Single Wall w/ Heating Base



► Key features:

- Rapid temperature change with heating base
- Temperature control up to 90 °C
- Designed for dominated strains of **microbes**



FS-V-D05L

▶ FS-V-SA05 Type Vessel-**Thermostat** Double Jacketed

▶ Key features:

- Special for low water level culture conditions
- Broken type, Anchor type and Spiral type impellers available.
- Solid state vessel performs 0-90 degrees rotation and 120 degrees for harvest.
- Ideal for **filamentous fungi**



Broken



Anchor




Spiral



5L solid state

► Customization for Solid State Cultures

➤ 3L Solid State Vessel



Front view

Rear view




#CUSTOMIZATION

Four kinds of impellers available!



► Customization for Solid State Cultures

➤ Humidifier for Solid State Cultures



#CUSTOMIZATION




Control the humidity
to get fruit body


Relative Humidity





► Specifications





	Vessel Type	Double Jacketed Dish Bottom Vessel (FS-V-A series)					
	Material	Borosilicate glass / 316L stainless steel for headplate and all fittings					
	Working Volume **	500ml	1L	3L	5L	10L	
	Total Volume Δ	1L	1.5L	3.8L	6.8L	12.5L	

	Vessel Type	Single Wall Dish Bottom Vessel (FS-V-B series)			
	Material	Borosilicate glass / 316L stainless steel for headplate and all fittings			
	Working Volume **	1L	3L	5L	10L
	Total Volume Δ	1.5L	3.8L	6.8L	12.5L

	Vessel Type	Air Lifter Vessel (FS-V-C series)			
	Material	Borosilicate glass / 316L stainless steel for headplate and all fittings			
	Working Volume **	5L single wall		5L double jacketed	
	Total Volume Δ	7L			

	Vessel Type	Single Wall Dish Bottom Vessel With Heating Blanket (FS-V-B series)					
	Material	Borosilicate glass / 316L stainless steel for headplate and all fittings					
	Working Volume **	1L	3L	5L	10L	15L	20L
	Total Volume Δ	1.5L	3.8L	6.8L	12.5L	18.7L	23.7L

	Vessel Type	Single Wall Plain Bottom Vessel With Heating Base Unit (FS-V-D series)				
	Material	Borosilicate glass / 316L stainless steel for headplate and all fittings				
	Working Volume **	3L		5L		10L
	Total Volume Δ	3.7L		6.7L		13.1L

	Vessel Type	Solid State (FS-V-SA05P)			
	Material	Borosilicate glass / 316L stainless steel for headplate and all fittings			
	Working Volume **	5L			
	Total Volume Δ	6.8L			

** Suggested Max.

Δ Total volumes are approximate and may vary slightly.

► Vessel Applications



Vessel Application

Application \ Vessel	FS-V-A series	FS-V-B series	FS-V-C series	FS-V-B series	FS-V-D series	FS-V-SA05P
	Double Jacketed Dish Bottom Vessel	Single Wall Dish Bottom Vessel	Air Lifter Vessel	Single Wall Dish Bottom Vessel with Heating Blanket	Single Wall Plain Bottom Vessel with Heating Base Unit	Solid State
Mammalian Cell Culture	● ●	● ○	○ ○	● ○	○ ○	○ ○
Aerobic Microorganism Culture (Note 1)	● ●	● ●	● ●	● ●	● ●	○ ○
Micro-aerobic Microorganism Culture (Note 2)	● ●	● ●	○ ○	● ●	● ●	○ ○
Anaerobic Microorganism Culture (Note 3)	● ●	● ●	○ ○	● ●	● ●	○ ○
Fragile Cell Culture (Note 4)	● ●	● ○	● ●	● ○	○ ○	○ ○
Photosynthesis Cell Culture (Note 5)	● ○	● ●	● ●	○ ○	● ○	○ ○
Plant Cell Culture	● ○	● ○	● ●	○ ○	○ ○	○ ○
Insect Cell Culture	● ●	● ○	○ ○	● ○	○ ○	○ ○
Solid State / Semi-solid State	○ ○	○ ○	○ ○	○ ○	○ ○	● ●

● ● Excellent

● ○ Good

○ ○ Not recommended

Note:

1. Some bacteria; yeast; fungi

2. Facultative culture (For example, some Lactobacillus; ethanol production)

3. Same as Note 2

4. This vessel is excellent for fragile cells that are easily sheared by any type of mechanical impeller


5. Plant; algae; cyanobacteria (blue-green algae)

Topics



- ▶ Benchtop fermentation system
 - Controller
 - Vessel Types
 - Optional Devices & Accessories
 - Fermenter Interface & Features


► Optional Devices



► Optional Devices

To enhance fermentation control experience

- Oxygen Enrichment Module (FS-O-OE)
 - Supply O₂ to Maintain aerobic environment
 - Flow rate duration via solenoid valve
 - Manually flow rate adjustment via flow meter




► Optional Devices

To enhance fermentation control experience

➤ Oxygen Enrichment Module with Mass Flow Control (FS-O-MF)

- Supply O₂ to Maintain aerobic environment
- More precise and accurate control
- Automatic control for gas flow



► Optional Devices

To enhance fermentation control experience

➤ Gas Mixing Station (FS-O-GM)

- Supply air, O₂, N₂ and CO₂ independently
- O₂ to increase DO level
- N₂ to create anaerobic environment
- CO₂ as alternative for pH control or provide carbonate for cell culture



► Optional Devices

To enhance fermentation control experience



- CO₂/O₂ Off-Gas Analyzer (FS-O-
 - Measure the amount of CO₂/O₂ off- gas
 - Monitor metabolic activities that exhale CO₂/O₂




**Methane Off-Gas Analyzer
(*Customized)**


► Optional Devices

To enhance fermentation control experience 




➤ Online Cell Density Devices

- Measure the amount Monitor cell growth rate without sampling
- Reduce the risk of contamination form sampling 
- Integration with fermentation system



Total Cell Density / Turbidity
Measurement at NIR



Viable Cell Density
Permittivity

► Optional Devices

To enhance fermentation control experience

➤ External Pump

- Digital control by microprocessor controller, more precise
- Easy load pump head
- Expand with the 2nd pump head with MU-D01/02



MU-D01/ MU-D02



MU-D03



Expand your feeding pumps!

► **Optional Devices**
To enhance fermentation control experience




- **Photo-bioreactor Lighting Module (FS-O-PB)**
 - Artificial light for photosynthesis
 - Adjustable light intensity
 - Fluorescent light source



► Application of Off-Gas Analyzer & Software Customization



➤ Get the information about cell metabolism via **CER** and **OUR**



✧ Monitor culture status to improve production process and yields



CER (Carbon Dioxide Evolution Rate) - the rate that CO₂ is being consumed
OUR (Oxygen Uptake Rate) - the rate that is O₂ being produced




Topics



► Benchtop fermentation system

- Introduction
- Vessel Types
- Optional Devices
- Fermenter Interface & Features

► Fermenter Interface & Features ► Intuitive user interface



► Fermenter Interface & Features

➤ Charting



► Optional Accessories



**Motor Shaft
Protection Cap**



**Fermentation
Bottle Holder**



Sampling Device



**Stainless Steel
Supporting foot**



**Composite Vessel
Handles**

Vessel Stand



Headplate Stand

Many Thanks!

Questions?