



### Precise and Intuitive Control

Every detail of the Integra 3 is geared towards one goal – to ensure you create the best possible conditions for ICSI and sample manipulation. When performing ART micromanipulation procedures, accuracy is paramount. Integra 3 offers the finest possible movement of the XY stage, precise movement of the micropipette through the fine and coarse control levers, and a choice of oil or air syringes for the control of specimens. Its further design features optimise ease of use.

- **Super-Fast** - Improved rapid micropipette set-up with one-step angle adjustment
- **Easy To Use** - Simple touch screen controls for adjusting stage temperature
- **Motion-Triggered Light** - A motion triggered light illuminates under the stage as soon as it detects your hand, to help you see when changing objectives
- **Shortcuts** - If you need to keep track of samples or procedures, just hit the shortcut buttons to film or photograph any stage of the process. There is also a counter button and an on-screen stopwatch
- **Safe** - Stage temperature control status indicator light and alarm monitors every second of the process. A log is kept of every alarm triggered
- **Colour Touch Screen** - 3.5 inch touch screen display showing temperature, stopwatch and counter. The electronic height indicators track the vertical position of the toolholders.
- **Integrated** - Large built-in mechanical XY stage plus temperature controllers with easy access buttons and touch screen