

RI Viewer™ Intuitive Software as Standard

In addition to its intuitive use and safety assurances, the new Saturn 5 features RI Viewer software. The software offers uncluttered full screen imaging from the microscope, digital magnification and a modern, clean user interface. RI Viewer also offers a recording function, line measurement tools which are visible on-screen and are printable, plus a built-in laser simulator for training and demonstration purposes. RI Viewer supports a range of languages.

- **Precise** - Sub-micron accuracy and unique computer controlled laser with guaranteed laser alignment
- **Rapid** - Faster than ever directional laser increases functionality and decreases procedure times
- **Always Spot-On** - As the Saturn 5 pilot laser travels down the same fibre optic path as the ablation laser, it guarantees consistent positioning
- **Foot Pedal** - A USB single or triple foot pedal can be programmed to control main features such as Fire, Zoom, Take photo, etc
- **CE Marked and FDA Approved** - In the USA, FDA cleared for clinical use for embryo biopsy, trophoctoderm biopsy, blastocyst collapsing and Laser Assisted Hatching (LAH)

